

## 7. Ordering Information

Components for the HPVA Series analyzers can be ordered using one of the following methods:

- Call Customer Service at (770) 662-3636
- Contact your local sales representative

When ordering, please use the information provided below to place your order.

Part Number	Item and Description
P02-25805-00	Gas inlet line, stainless steel, 8 ft.
P02-31701-00	Furnace insulator top, half
003-51123-00	Fuse, 6.25 Amp, 3AG, Slow Blow
003-51134-01	Fuse, 5.0 Amp, 5x20 mm, Slow Blow
004-25040-01	Gasket, 1/4 in., stainless steel, silver plated; for sample cylinder
004-27056-01	Filtered gasket, 60 micron; for sample cylinder when fine powders used
P02-25821-00	Sample cylinder, 2 cc
P02-25827-00	Sample cylinder, 10 cc
P02-25846-00	Sample tube funnel
004-25466-02	O-ring -010 60 Duro Viton F, for degas or analysis port
P02-25820-00	Sample isolation valve
P02-25823-00	Sample tube stem, used with sample isolation valve
P02-33012-00	Cryogenic Kit, includes Dewar and accessories
P02-34013-00	Sample Chamber Assembly, 0.75cc
004/16847/00	5A Zeolite High Pressure Reference Material Kit