

Consumables

Consumables are listed below. Contact your Particulate Systems representative to reorder consumables according to the exchange limit and life span information shown.

| Item | Approximates Exchange Limit or Life Span |
|----------------------|--|
| Tubes (1) and (2) | Every exchange of titration solution |
| Tubes (3), (5) - (9) | Every exchange of sample solution |
| Tube (4) | 1 month |
| Syringe packing | 6 months |

| Part Number | Item and Description |
|-----------------------|--------------------------------|
| P09-00003-00002 | pH electrode, AutoTitrator |
| P09-00003-00003 | Stirrer, pH; AutoTitrator |
| Size Flow Cell | |
| P09-00000-04001 | Quartz Cuvette, size flow cell |
| P09-00000-04002 | Tubing, Teflon, size flow cell |
| P09-00000-04003 | Luer fitting, size flow cell |