

SPECIFICATIONS

Particle Insight



PHYSICAL

Height	38.1 cm (15 in.)
Width	25.4 cm (10 in.)
Depth	63.5 cm (25 in.)
Weight	13.2 kg (29lbs.)

ELECTRICAL

Voltage	100 - 240 VAC.
Frequency	50 - 60 Hz.

ENVIRONMENT

Temperature	10 to 45 °C (50 to 113 °F.) -10 to 55 °C (14 to 131 °F.)
-------------	---

COMPUTER REQUIREMENTS

Processor	2 GHz CPU or superior
Memory	1 gigabyte of RAM.
Hard Dish Space	60 gigabytes.
Monitor	58.42 cm (23 in.) Flat Panel.
Media Drive	Writable CD or DVD.
Interface	IEEE 1394 (Firewire "A" type.)
Software	Windows® 7. Windows XP.

Due to continuous improvements, specifications are subject to change without notice.