

SPECIFICATIONS

FORMULATION - TURBISCAN AGS



TECHNOLOGY

S-MLS 880 NM

Sample Cell Volume	20 mL
Temperature range	RT - 60°C (3 racks) *90°C on demand
Number of Samples	54
Maximum sample conc.	95% v/v
Measured size range	10 nm - 1mm
Reproducibility/Repeatability	0.1 % / 0.05%
Acquisition scan step	40µm
Automatic samples recognition (bar-code)	Yes
ISO Certification	TR 13097 Compliant
Dimensions	145 x 75 x 85 cm
Weight	50kg