

# SPECIFICATIONS

## FORMULATION - TURBISCAN CLASSIC 2



### TECHNOLOGY

### S-MLS 850 NM

Sample Cell Volume	7mL
Temperature range	RT
Number of Samples	1
Maximum sample conc.	60% v/v
Measured size range	10 nm - 1mm
Reproducibility/Repeatability	0.25% / 0.1%
Acquisition scan step	20 $\mu$ m
Automatic samples recognition (bar-code)	No
ISO Certification	TR 13097 Compliant
Dimensions	34 x 16 x 29 cm
Weight	5kg