

# SPECIFICATIONS

## FORMULATION - TURBISCAN LAB



### TECHNOLOGY

### S-MLS 880 NM

Sample Cell Volume	4 or 20 ml
Temperature range	RT - 60°C
Number of Samples	1
Maximum sample conc.	95% v/v
Measured size range	10 nm - 1mm
Reproducibility/Repeatability	0.1% / 0.05%
Acquisition scan step	40µm
Automatic samples recognition (bar-code)	Yes
ISO Certification	TR 13097 Compliant
Dimensions	38 x 42 x 32 cm
Weight	13kg