

SPECIFICATIONS

FORMULATION - TURBISCAN TOWER



TECHNOLOGY

S-MLS 880 NM

Sample Cell Volume	4 or 20 mL
Temperature range	4 - 80°C
Number of Samples	6
Maximum sample conc.	95% v/v
Measured size range	10 nm - 1mm
Reproducibility/Repeatability	0.05% / 0.05%
Acquisition scan step	20µm
Automatic samples recognition (bar-code)	Yes
ISO Certification	TR 13097 Compliant
Dimensions	35 x 45 x 90 cm
Weight	45kg